



Anti-FKBP15 (N-terminal) polyclonal antibody (CABT-BL1518)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to N-terminus 15 amino acids of human FKBP15.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Format	Liquid
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Entrez Gene ID	23307
Protein Refseq	NP_056073
UniProt ID	Q5T1M5