



Anti-FKH2 (aa 611-710) polyclonal antibody (CABT-BL1519)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: SDSADKSTNNNGGKMNLP AISTSSLDENG NLEPTTTTSS GNSNSVPQTGTTTSSLAANSLRLSQPYDTLLRSPTKAFHI TAMEAYTPERGSANRARSPL, corresponding to amino acids 611/710 of <i>S. cerevisiae</i> Fkh2p.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>Saccharomyces cerevisiae</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [855656](#)

Protein Refseq [NP_014331](#)

UniProt ID [P41813](#)
