



Anti-Fla2 (internal region) polyclonal antibody (CABT-BL1520)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from an internal domain of Flagellin 2 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Clostridia Flagellum
Conjugate	Unconjugated
Applications	WB
Format	Antiserum
Preservative	None
Storage	Short term storage: +4°C. Long term storage: -20°C.

BACKGROUND

Introduction	AF333971 Arabidopsis thaliana fasciclin-like arabinogalactan-protein 2 (Fla2) mRNA, complete cds
--------------	--

GENE INFORMATION

Entrez Gene ID	3728773
Protein Refseq	YP_359843

UniProt ID

[Q3ADE1](#)
