



Anti-FMO1 (aa 1-407) polyclonal antibody (CABT-BL1532)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 407 of Human FMO1 (Q01740)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Microsome membrane. Endoplasmic reticulum membrane.
Format	Liquid
Buffer	10% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	catalyzes thiobenzamide S-oxidation; has additional S- and N-oxidation activities as well
---------------------	---

GENE INFORMATION

Entrez Gene ID [2326](#)

Protein Refseq [NP_002012](#)

UniProt ID [Q01740](#)
