



Anti-FNDC1 polyclonal antibody (CABT-BL1534)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human FNDC1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	May be an activator of G protein signaling.
--------------	---

GENE INFORMATION

Entrez Gene ID [84624](#)

Protein Refseq [NP_115921.2](#)

UniProt ID [Q4ZHG4](#)
