



Anti-FNDC3A (full length) polyclonal antibody (CABT-BL1536)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length FNDC3A protein (Human), (ENSP00000347952, 1134 a.a).
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Fibronectin type-III domain-containing protein 3a is a protein that in humans is encoded by the FNDC3A gene.
--------------	--

GENE INFORMATION

Entrez Gene ID [22862](#)

Protein Refseq [NP_055738](#)

UniProt ID [A0A024RDT9](#)
