



Anti-FOLH1B (aa 160-190) polyclonal antibody (CABT-BL1539)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 160-190 of Human FOLH1B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Cellular Localization	Cytoplasmic
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	219595
Protein Refseq	NP_710163

UniProt ID

[Q9HBA9](#)
