



Anti-FOXD3 polyclonal antibody (CABT-BL1544)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from internal of Human FOXD3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nucleus.
Format	Liquid
Buffer	50% Glycerol, PBS
Preservative	0.03% Sodium Azide
Storage	Store at 4°C. Aliquot and store at -20°C long-term.

BACKGROUND

Introduction	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. Mutations in this gene cause autoimmune susceptibility 1.
---------------------	---

GENE INFORMATION

Entrez Gene ID [27022](#)

Protein Refseq [NP_036315](#)

UniProt ID [Q9UJU5](#)
