



Anti-Chicken FOXD3 polyclonal antibody (CABT-BL1545)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from N terminal domain of chicken FOXD3 protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Chicken
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	IHC, WB
Cellular Localization	Nucleus.
Format	Liquid
Size	100 µl
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. Mutations in this gene cause autoimmune susceptibility 1.
---------------------	---

GENE INFORMATION

Entrez Gene ID	395794
Protein Refseq	NP_990282
UniProt ID	P79772