



# Anti-FoxD3 polyclonal antibody (CABT-BL1547)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB, ICC
Preservative	None
Storage	-20°C, Avoid freeze/thaw

## BACKGROUND

Introduction	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. Mutations in this gene cause autoimmune susceptibility 1.
--------------	---

## GENE INFORMATION

Entrez Gene ID	<a href="#">15221</a>
Protein Refseq	<a href="#">NP_034555</a>
UniProt ID	<a href="#">A2BDY3</a>