



Anti-FOXD3 polyclonal antibody (CABT-BL1548)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB, ICC
Preservative	None
Storage	2-8°C (up to 1 month) below -15°C (long term); Avoid repeat freeze thaw cycles

BACKGROUND

Introduction	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. Mutations in this gene cause autoimmune susceptibility 1.
--------------	---

GENE INFORMATION

Entrez Gene ID	27022
Protein Refseq	NP_036315
UniProt ID	Q9UJU5