



Anti-FOXF2 polyclonal antibody (CABT-BL1550)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB
Preservative	None
Storage	2-8°C (up to 1 month) below -15°C (long term); Avoid repeat freeze thaw cycles

BACKGROUND

Introduction	FOXF2 encodes forkhead box F2, one of many human homologues of the Drosophila melanogaster transcription factor forkhead. FOXF2 is expressed in lung and placenta, and has been shown to transcriptionally activate several lung-specific genes.
--------------	--

GENE INFORMATION

Entrez Gene ID	2295
Protein Refseq	NP_001443
UniProt ID	Q12947