



Anti-FOXI3 (aa 247-260) polyclonal antibody (CABT-BL1554)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: RLRVSGKLEGDSPS by a Cysteine residue linker, corresponding to internal sequence amino acids 247-260 of Mouse FOXI3.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse, Rat
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Size	200 µl
Buffer	99% Tris buffered saline, 0.5% BSA pH: 7.30
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction FOXI3 (forkhead box I3) is a protein-coding gene. GO annotations related to this gene include

double-stranded DNA binding and sequence-specific DNA binding.

GENE INFORMATION

Entrez Gene ID	232077
----------------	------------------------

Protein Refseq	NP_001094934
----------------	------------------------------

UniProt ID	D3Z120
------------	------------------------
