



Anti-FOXI3 (aa 396-445) polyclonal antibody (CABT-BL1555)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 396 - 445 (EASADTLQLS NSTSNSTGQR SSYYSPFPAS TSGGQSSPFS SPFHNFSMVN) of human FOXI3 (XP_001126001).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	FOXI3 (forkhead box I3) is a protein-coding gene. GO annotations related to this gene include double-stranded DNA binding and sequence-specific DNA binding.
---------------------	--

GENE INFORMATION

Entrez Gene ID [344167](#)

Protein Refseq [NP_001129121](#)

UniProt ID [A8MTJ6](#)
