



Anti-FOXL2 (aa 364-376) polyclonal antibody (CABT-BL1557)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: DSKTGALHSRLDL, corresponding to amino acids 364-376 of FOX L2.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Cellular Localization	Nucleus.
Format	Liquid
Size	100 µg
Buffer	0.5% BSA, Tris-saline. pH 7.3
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction This gene encodes a forkhead transcription factor. The protein contains a fork-head DNA-

binding domain and may play a role in ovarian development and function. Mutations in this gene are a cause of blepharophimosis syndrome and premature ovarian failure 3.

GENE INFORMATION

Entrez Gene ID	668
----------------	---------------------

Protein Refseq	NP_075555
----------------	---------------------------

UniProt ID	P58012
------------	------------------------
