



Anti-FOXL2 polyclonal antibody (CABT-BL1559)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide surrounding amino acid 367 of human FOXL2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	WB
Preservative	None
Storage	-20°C

BACKGROUND

Introduction	This gene encodes a forkhead transcription factor. The protein contains a fork-head DNA-binding domain and may play a role in ovarian development and function. Mutations in this gene are a cause of blepharophimosis syndrome and premature ovarian failure 3.
--------------	--

GENE INFORMATION

Entrez Gene ID	668
----------------	---------------------

Protein Refseq [NP_075555](#)

UniProt ID [P58012](#)
