



Anti-FTHL17 (aa 38-87) polyclonal antibody (CABT-BL1571)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region within N-terminal amino acids 38-87 (MAFYFNRRDV ALENFFRYFL RLSDDKMEHA QKLMRLQNLR GGHICLHDIR) of Human FTHL17.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	This gene is orthologous to a mouse gene that encodes a ferritin heavy polypeptide-like protein and shows testis-specific expression.
---------------------	---

GENE INFORMATION

Entrez Gene ID [53940](#)

Protein Refseq [NP_114100](#)

UniProt ID [Q9BXU8](#)
