



Anti-FUZ (aa 140-224) polyclonal antibody (CABT-BL1573)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human FUZ.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Probable planar cell polarity effector involved in cilium biogenesis. May regulate protein and membrane transport to the cilium. May regulate the morphogenesis of hair follicles which depends on functional primary cilia.
--------------	--

GENE INFORMATION

Entrez Gene ID [80199](#)

Protein Refseq [NP_001165408](#)

UniProt ID [Q9BT04](#)
