



# Anti-G6PC2 (aa 150-250) polyclonal antibody (CABT-BL1577)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide conjugated to KLH derived from within residues 150 - 250 of Mouse IGRP.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Cell Membrane
<b>Format</b>	Liquid
<b>Buffer</b>	1% BSA, PBS, pH 7.4
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	Glucose-6-phosphatase 2 is an enzyme that in humans is encoded by the G6PC2 gene.
---------------------	---

# GENE INFORMATION

Entrez Gene ID [14378](#)

---

Protein Refseq [NP\\_067306](#)

---

UniProt ID [Q9Z186](#)

---