



Anti-GABARAPL1 polyclonal antibody (CABT-BL1578)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Human GABARAPL1 recombinant protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasm. Membrane. Endoplasmic reticulum. Golgi apparatus.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Ubiquitin-like modifier that increases cell-surface expression of kappa-type opioid receptor through facilitating anterograde intracellular trafficking of the receptor. Involved in formation of autophagosomal vacuoles. Whereas LC3s are involved in elongation of the phagophore membrane, the GABARAP/GATE-16 subfamily is essential for a later stage in autophagosome
---------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

maturation.

GENE INFORMATION

Entrez Gene ID	23710
----------------	-----------------------

Protein Refseq	NP_113600
----------------	---------------------------

UniProt ID	A0A024RAP5
------------	----------------------------
