



Anti-GAK polyclonal antibody (CABT-BL1589)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to N-terminus of human GAK.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein G purification
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Format	Liquid
Buffer	In PBS (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Associates with cyclin G and CDK5. Seems to act as an auxilin homolog that is involved in the uncoating of clathrin-coated vesicles by Hsc70 in non-neuronal cells. Expression oscillates slightly during the cell cycle, peaking at G1.
---------------------	--

GENE INFORMATION

Entrez Gene ID [2580](#)

Protein Refseq [NP_005246](#)

UniProt ID [B4DS37](#)
