



Anti-GALNTL4 (internal region) polyclonal antibody (CABT-BL1595)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide from a region within internal amino acids 468-517 (IAYGVLQNSL KTDLC LDQGP DTENV PIMYI CHGMTPQNVY YTSSQ QIHVG) of Human GALNTL4
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Golgi Apparatus membrane; Single-pass type II membrane protein
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	374378
Protein Refseq	NP_940918

UniProt ID

[Q58A54](#)
