



Anti-GATS (full length) polyclonal antibody (CABT-BL1597)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant full length human STAG3OS (AAH01129, 1 a.a. ~ 102 a.a).
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Buffer	1X PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	352954
Protein Refseq	NP_849153
UniProt ID	Q8NAP1