



Anti-GBP7 (aa 530-560) polyclonal antibody (CABT-BL1598)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within of C terminal amino acids 530-560 of Human GBP7.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane; Multi-pass membrane protein
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	Guanylate-binding proteins, such as GBP7, are induced by interferon and hydrolyze GTP to both GDP and GMP (Olszewski et;al., 2006 (PubMed 16689661)).
---------------------	---

GENE INFORMATION

Entrez Gene ID [388646](#)

Protein Refseq [NP_997281](#)

UniProt ID [Q8N8V2](#)
