



Anti-gcy-23 (aa 821-920) polyclonal antibody (CABT-BL1608)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: MDYIRDGKTKFRPTVHDRTQIHPDLVALLDCWNENPEVR PSIRRVRLNTENYLKVKGSLVDQMMRMMEQYANNLEKLVA ERTGMLEEANVRADKLLGQL, corresponding to amino acids 821/920 of <i>Caenorhabditis elegans</i> Gcy23
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>Caenorhabditis elegans</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	191652
Protein Refseq	NP_500309

UniProt ID

[G5EEE9](#)
