



Anti-GEN1 (aa 67-97) polyclonal antibody (CABT-BL1618)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH conjugated synthetic peptide corresponding to a region within N terminal amino acids 67-97 of Human GEN1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	348654
Protein Refseq	NP_001123481

UniProt ID

[Q17RS7](#)
