



Anti-GET4 (aa 277-327) polyclonal antibody (CABT-BL1619)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------------|--|
| Immunogen | Synthetic peptide, corresponding to a region within amino acids 277-327 of Human CEE |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | Immunogen affinity purified |
| Conjugate | Unconjugated |
| Applications | WB, IP |
| Format | Liquid |
| Size | 100 µl |
| Buffer | 0.1% BSA, Tris buffered saline |
| Preservative | 0.09% Sodium Azide |
| Storage | Store at +4°C. |

BACKGROUND

Introduction UPF0363 protein C7orf20 is a protein that in humans is encoded by the C7orf20 gene.

GENE INFORMATION

Entrez Gene ID [51608](#)

Protein Refseq [NP_057033](#)

UniProt ID [Q7L5D6](#)
