



Anti-GL (N-terminal) polyclonal antibody (CABT-BL1628)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from the N terminal domain of DCTN1 (Drosophila melanogaster)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr, WB
Cellular Localization	Cytoplasm. Cytoplasm; cytoskeleton.
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	39536
Protein Refseq	NP_524061

UniProt ID

[P13496](#)
