



Anti-GLYATL2 (aa 109-158) polyclonal antibody (CABT-BL1634)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 109-158 (EGLDEAIRKV ATSKSVQVDY MKTILFIPEL PKKHKTSSND KMELFEVDDD) of Human GLYL2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Mitochondrial
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Entrez Gene ID [219970](#)

Protein Refseq [NP_659453](#)

UniProt ID [A0A024R4Z5](#)
