



Anti-GM5581 (aa 20-69) polyclonal antibody (CABT-BL1636)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 20-69 (EEWECLDSAQ RDLYRDVMLE NYHNLVSVG AVSKSEVIFC LEQNNESWIA) of Mouse Gm5581 (XP_485832).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	98% PBS, 2% Sucrose
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [434094](#)