



Anti-GM6185 (aa 700-749) polyclonal antibody (CABT-BL1637)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 700-749 (MYSITKGYVK SQEDAKLLIR QISGRESIYR KLYDILEKNK REAIRELGLI) of Mouse Gm6185 (XP_901143.2).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [620839](#)