



# Anti-GOLGA6B polyclonal antibody (CABT-BL1641)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugate	Unconjugated
Applications	WB, ELISA
Buffer	PBS with 0.09% (W/V) sodium azide
Preservative	0.09% Sodium Azide
Storage	GOLGA6B Antibody (Center) can be refrigerated at 2-8°C for up to 6 months. For long term storage, place the at -2°C.

## GENE INFORMATION

Entrez Gene ID	<a href="#">55889</a>
Protein Refseq	<a href="#">NP_061122.4</a>
UniProt ID	A6NDN3