



# Anti-GOSR1 (aa 2-155) polyclonal antibody (CABT-BL1646)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |  |
|--------------------|--|
| Immunogen          | Recombinant cytoplasmic domain of rat GOSR2(aa 2 - 155). |
| Isotype            | IgG  |
| Source/Host        | Rabbit   |
| Species Reactivity | Rat  |
| Conjugate          | Unconjugated   |
| Applications       | WB, ICC  |
| Preservative       | None   |
| Storage            | -20°C, Avoid freeze / thaw cycles                        |

## BACKGROUND

|              |  |
|--------------|--|
| Introduction | This gene encodes a trafficking membrane protein which transports proteins among the endoplasmic reticulum and the Golgi and between Golgi compartments. This protein is considered an essential component of the Golgi SNAP receptor (SNARE) complex. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. |
|--------------|--|

## GENE INFORMATION

|                |                       |
|----------------|-----------------------|
| Entrez Gene ID | <a href="#">94189</a> |
|----------------|-----------------------|

Protein Refseq [NP\\_446036](#)

---

UniProt ID [Q62931](#)

---