



# Anti-GOT1L1 (aa 353-383) polyclonal antibody (CABT-BL1648)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide conjugated to KLH, corresponding to a region within the C terminal amino acids 353-383 of Human GOT1L1.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Cytoplasm.
<b>Format</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

Entrez Gene ID [137362](#)

Protein Refseq [NP\\_689626](#)

---

UniProt ID [Q8NHS2](#)

---