



Anti-GPR101 (C-terminal) polyclonal antibody (CABT-BL1654)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthesized peptide derived from C-terminal region of human GPCR GPR101.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, ICC/IF
Cellular Localization	Cell membrane.
Format	Liquid
Buffer	50% Glycerol, PBS, 150mM Sodium chloride, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Orphan receptor.
---------------------	------------------

GENE INFORMATION

Entrez Gene ID [83550](#)

Protein Refseq [NP_473362](#)

UniProt ID [Q96P66](#)
