



# Anti-GPR113 polyclonal antibody (CABT-BL1656)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human GPR113.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Probable G-protein coupled receptor 113 is a protein that in humans is encoded by the GPR113 gene.
--------------	--

## GENE INFORMATION

Entrez Gene ID [165082](#)

---

Protein Refseq [NP\\_001138640](#)

---

UniProt ID [Q8IZF5](#)

---