



Anti-GPR133 polyclonal antibody (CABT-BL1657)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide (Human) conjugated to KLH
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Cellular Localization	Cell membrane; Multi-pass membrane protein.
Format	Liquid
Buffer	PBS
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Orphan receptor.
--------------	------------------

GENE INFORMATION

Entrez Gene ID [283383](#)

Protein Refseq [NP_942122](#)

UniProt ID [Q6QNK2](#)
