



Anti-GPR133 polyclonal antibody (CABT-BL1658)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to human GPR133.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	Immunoaffinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Orphan receptor.
--------------	------------------

GENE INFORMATION

Entrez Gene ID [283383](#)

Protein Refseq [NP_942122](#)

UniProt ID [Q6QNK2](#)
