



# Anti-GPR172B (internal region) polyclonal antibody (CABT-BL1661)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from an internal domain of mouse GPR172B protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell Membrane; Multi-pass membrane protein
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

Introduction	GPR172B acts as receptor for porcine endogenous retrovirus subgroup A and is strongly expressed in testis.
--------------	--

# GENE INFORMATION

Entrez Gene ID [52710](#)

---

Protein Refseq [NP\\_083919](#)

---

UniProt ID [Q9D8F3](#)

---