



Anti-GPR172B (full length) polyclonal antibody (CABT-BL1662)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length human protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein G purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell Membrane; Multi-pass membrane protein
Format	Liquid
Buffer	1X PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	GPR172B acts as receptor for porcine endogenous retrovirus subgroup A and is strongly expressed in testis.
---------------------	--

GENE INFORMATION

Entrez Gene ID [55065](#)

Protein Refseq [NP_001098047](#)

UniProt ID [Q9NWF4](#)
