



# Anti-GPR6 (aa 1-50) polyclonal antibody (CABT-BL1663)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide corresponding to a region within N terminal amino acids 1-50 (NASAASLNDS QVVVVAAEGA AAAATAAGGP DTGEWGPPAA AALGAGGGAN) of Human GPCR GPR6.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ICC/IF
<b>Cellular Localization</b>	Cell membrane; Multi-pass membrane protein.
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	2% Sucrose, PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction** G protein-coupled receptor 6, also known as GPR6, is a protein which in humans is encoded by

the GPR6 gene.

---

## GENE INFORMATION

---

Entrez Gene ID	<a href="#">2830</a>
----------------	----------------------

---

Protein Refseq	<a href="#">NP_005275</a>
----------------	---------------------------

---

UniProt ID	<a href="#">P46095</a>
------------	------------------------

---