



Anti-GRAMD1A (aa 266-296) polyclonal antibody (CABT-BL1677)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 266-296 of Human GRAM1DA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	99% PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	57655
Protein Refseq	NP_065946

UniProt ID

[Q96CP6](#)
