



Anti-groES (full length) polyclonal antibody (CABT-BL1686)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length GroES purified chromatographically from an E.coli strain overexpressing the protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	E. coli
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IP
Cellular Localization	Cytoplasmic
Format	Liquid
Size	50 µl
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [959981](#)

Protein Refseq [NP_290775](#)

UniProt ID [A0A075UCT7](#)
