



Anti-GSDMB (internal region) polyclonal antibody (CABT-BL1690)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------------|--|
| Immunogen | This antibody is specific for the Middle Region of the target protein. |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | Affinity purified |
| Conjugate | Unconjugated |
| Applications | WB, ELISA, IHC-P |
| Buffer | 20 mM Potassium Phosphate, 150 mM Sodium Chloride, pH 7.0 |
| Preservative | None |
| Storage | Store at -20°C. Avoid freeze-thaw cycles. |

BACKGROUND

| | |
|--------------|--|
| Introduction | This gene encodes a member of the gasdermin-domain containing protein family. Other gasdermin-family genes are implicated in the regulation of apoptosis in epithelial cells, and are linked to cancer. Multiple transcript variants encoding different isoforms have been found for this gene. Additional variants have been described, but they are candidates for nonsense-mediated mRNA decay (NMD) and are unlikely to be protein-coding. |
|--------------|--|

GENE INFORMATION

Entrez Gene ID [55876](#)

Protein Refseq [NP_001035936](#)

UniProt ID [Q8TAX9](#)
