



Anti-GTF3C1 (aa 1475-1525) polyclonal antibody (CABT-BL1705)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region between residue 1475 and 1525 of human TFIIC / GTF3C1 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IP
Cellular Localization	Nuclear
Format	Liquid
Buffer	8mM PBS, 60mM Citrate, 150mM Tris, pH 7-8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	The transcription factor IIIC (TFIIIC) is a multi-subunit DNA-binding factor required for promoter recognition and TFIIIB assembly on tRNA genes transcribed by RNA polymerase III.
---------------------	---

GENE INFORMATION

Entrez Gene ID [2975](#)

Protein Refseq [NP_001511](#)

UniProt ID [Q12789](#)
