



Anti-G γ -2-Subunit polyclonal antibody (CABT-BL1716)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	recombinant, human, retina Gy-2-subunit
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ELISA, IB
Format	Liquid.
Preservative	None
Storage	-20°C, Avoid freeze/thaw
Ship	Blue Ice Only Standard Handling

GENE INFORMATION

Protein Refseq	NP_444292
----------------	---------------------------