



Anti-HADHB polyclonal antibody (CABT-BL1748)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 321-351 of Human HADHB, conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Mitochondrion.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	This gene encodes the beta subunit of the mitochondrial trifunctional protein, which catalyzes the last three steps of mitochondrial beta-oxidation of long chain fatty acids. The mitochondrial membrane-bound heterocomplex is composed of four alpha and four beta subunits, with the
---------------------	--

beta subunit catalyzing the 3-ketoacyl-CoA thiolase activity. The encoded protein can also bind RNA and decreases the stability of some mRNAs. The genes of the alpha and beta subunits of the mitochondrial trifunctional protein are located adjacent to each other in the human genome in a head-to-head orientation. Mutations in this gene result in trifunctional protein deficiency. Alternatively spliced transcript variants encoding different isoforms have been described.

GENE INFORMATION

Entrez Gene ID	3032
----------------	----------------------

Protein Refseq	NP_000174
----------------	---------------------------

UniProt ID	P55084
------------	------------------------
