



Anti-HAPLN4 (aa 251-300) polyclonal antibody (CABT-BL1750)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 251-300 (NYGYRHNAEE RYDAFCFTSN LPGRVFFLKP LRPVPFSGAA RACAARGAAV) of Human HAPLN4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Secreted; extracellular space; extracellular matrix.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [404037](#)

Protein Refseq [NP_075378](#)

UniProt ID [Q86UW8](#)
