



Anti-HAPLN4 (aa 30-135) polyclonal antibody (CABT-BL1751)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment with tag: QRGRKKVVHV LEGESGSVVV QTAPGQVVSH RGGTIVLPCR YHYEAAAHGH DGVRLKWKV VDPLAFTDVF VALGPQHRAF GSYRGRAELQ GDGPGDASLV LRNVT, corresponding to amino acids 30-135 of Human HAPLN4
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Secreted; extracellular space; extracellular matrix.
Format	Liquid
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [404037](#)

Protein Refseq [NP_075378](#)

UniProt ID [Q86UW8](#)
