



Anti-HAS polyclonal antibody (CABT-BL1754)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from Cter domain of Hyaluronan synthase from <i>Paramecium bursaria</i> chlorella virus-1 (PBCV-1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Virus
Conjugate	Unconjugated
Applications	WB
Format	Antiserum
Preservative	None
Storage	Short term storage: +4°C. Long term storage: -20°C.

BACKGROUND

Introduction	Hyaluronan synthases (HAS) are membrane-bound enzymes which use UDP-alpha-N-acetyl-D-glucosamine and UDP-alpha-D-glucuronate as substrates to produce the glycosaminoglycan hyaluronan at the cell surface and extrude it through the membrane into the extracellular space.
---------------------	--

GENE INFORMATION

Entrez Gene ID	918007
Protein Refseq	NP_048446

UniProt ID

[Q84419](#)
